Pelvik: it consists of an advanced system of three vaginal cones involving different weight measures (white 30g, light pink 47g, bright pink 65g) in order to gain strength to the pelvic floor muscles in a practical, effective and autonomous approach, resembling the technique of the first expert in this field, Dr Kegel. The pelvic floor is the muscular and ligamentous complex responsible to preserve important organs such as bladder, uterus, intestine and the correct functionality of the urethra, vagina and of the rectus.

The impairment of the pelvic floor could be due to:
 -pregnancy
 -menopause
 -physical sport activities
 -obesity
 -chronic constipation
 -chronic cough
 -heavy physical exertion
 -possible results of previous surgeries
 -familial predisposition

These conditions could lead to:
 -urinary incontinence (from exertion)
 -vaginal prolapse (vesicles, uterus, rectus)
 -fecal incontinence
 -unsatisfactory and/or painful sexual intercourse

Thank to the increase of female longevity, the woman has higher chances to live a significant period of her life with serious pelvic floor dysfunctions that interferes on particular levels:
 -Personal (with her intimacy, ex: the use of diapers for incontinence or pessary for the prolapse followed by other issues that derive from it)
 -Social (difficulties to be away from home and maintain interpersonal relationships)
 -Sexual (personal and/or partner's displeasure).

Limiting personal physical and partner's wellbeing.

Vaginal cones Pelvik.
 These vaginal cones have been developed by specialists to help women to follow correct exercises of the pelvic

floor. The Pelvik vaginal cones consisting in the inside of an iron sphere which by moving emits light vibrations, after inserting them in the vagina creates a feedback sensation on the walls helping the woman to gain perception of the cone and more stimulation on the pelvic musculature facilitating the work and increasing the efficacy. The cone itself with its weight has the tendency to escape from the vagina due to the simple gravity law inducing the stimulus to contract the perivaginal muscles denying the exit and by combining vibrations and weight lead to an "automatic training".

The muscular contraction exercise and the progressive and gradual switch to heavier cones, gain efficient strength to the perineal musculature, and allow to evolve in a bright improvement of the pelvic tone in just few weeks denoting its benefits.

How to use it

Understanding its own pelvic floor by inserting a finger, washed and properly lubricated, in the vagina for at least 2 cm. To feel the vaginal walls tightening when it will be induced contraction to the pelvic musculature, similar to the stopping reflex of urination during the micturition. At the beginning it is probably difficult to know its pelvic muscles and to experience the vaginal contraction. Thank to Pelvik cones and to its vibrating characteristic it will be easier to learn recognize and use the musculature even in the first phases of the exercise.

Cleaning of the cone and introduction in the vagina:

IMPORTANT: accurately clean the hands and disinfect the vaginal Pelvik cone with specific agents or disinfectant soaps (ex: purell or clorexidine), washing it with water and drying it with a soft tissue, do no use brushes or sponges and try not to scratch the cone, before and after the use. With the usage of a lubricant based up on water, the cone has to be inserted 3 cm over the inside of the vagina, leaving the tail of the cone on the outside helping the possible removal.

First approach:

Begin the procedure with the lightest cone (white) and look for the best position preferred. The first times it is suggested to stress the pelvic musculature while lying on a supine position and subsequently changing posture standing up inducing the cone, for the gravitational force, to easily bypass due to a weakened pelvic floor.

Practicing with Pelvik

Following the original method of Dr. Kegel, it is important to learn how to reduce and raise the pelvic floor along with its musculature. Each phase we will perform has to be comfortable, even if it could result in a strenuous dealing with the muscular training.

We begin the exercise, after introducing the vaginal cone, the pelvic musculature contracts for around 2 seconds holding the vaginal cone. Progressing with the exercise routine, the spasms lasting 2 seconds will promptly advance up to over ten seconds:

Begin with few seconds of contractions and reoccur for 10 times, with breaks of 5/10 seconds from each tightening. Initially increase the intervals to the personal need; it is important to the contraction period to last equally during each set.

WARNING: do not abuse of the abdominal muscles or the gluteal muscles to hold the cone, instead focus on the pelvic musculature. In case of loosing control of the pelvic musculature refer to a specialist.

Note that during the workout it is possible to check the functionality of the own pelvic muscles, by inserting a finger during the constriction and perceiving the recovery of the cone and by observing through a mirror the tail of the vaginal cone Pelvik that shortens on each spasm. The training is recommended at least three times a week for around twelve weeks or after three weeks at the earliest from the disappearance of the symptoms, and a salutary use for the maintenance of the pelvic wealth achieved. The use of the vaginal cones for preventive purposes is advised to every women.

Advanced exercises:

The more time is taken care of is the pelvic musculature, the more time the contraction will last, and less breaks will be necessary to end the set of training.

Afterwards it is possible to switch to the next weight, Trying to follow some specific exercises with the cone in the vagina taking into consideration its own level while:

- walking for at least one minute;
- walking up and down the stairs;
- Jumping on the spot for at least thirty seconds;
- Coughing for ten to fifteen times;
- Groveling;
- Washing hands under cold water for one minute.

Overcoming these exercises, we can go through to the next cone of higher weight initially reducing the time of contraction and prolonging the breaks.

The periods of time to switch the level of the vaginal cones depends on each woman.

Fulfilled a good level of training, the ovum can be applied for a longer period of time throughout the day and even when following other activities, as an "automatic workout".

The continuous use is not recommended for over a consecutive hour.

WARNINGS:

Keep out of reach of children
 After twenty weeks of daily and correct treat of the vaginal cones with no evidence of symptom improvement, it is suggested to lay off the use and consult a specialist with the purpose to establish variations to the therapy.



The vaginal cones DO NOT have to be used:

- during on course pregnancy
- Before the sixth/eighth week postpartum
- During the menstrual cycle
- Before the complete healing of the vaginal epithelium of any origin
- After a surgical intervention
- During vaginal mucous infections or difficulties
- If scraped or injured.


The cones weight could vary of few grams from the weight subscribed on the box

In case of any health problem, refer to medical counseling. This scripture does not want and must not be replaced to the sanitary expert opinion. Before undertake any therapy we suggest to always consult a specialist.

Pelvik vaginal cones Diameter: 29mm. Medical silicone.

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